

Core Imaging Facility

Ben Fowler, Manager (405) 271-7545 http://Imaginglinux2.omrf.ouhsc.edu/

Oklahoma Medical Research Foundation

825 N.E. 13th Street Oklahoma City, OK 73104 405-271-6673 800-522-0211 http://OMRF.org/

Brief Immunofluorescence Protocol

(Non-conjugated primary and conjugated secondary)

- 1. Remove embedding matrix
 - a) paraffin run through deparaffinization process
 - b) cryo rinse in water
- 2. Dry slightly (don't let tissue dry out; can be done with a KimWipe) and <u>circle tissue with PapPen</u>. Do this for one slide at a time.
- 3. Add PBS to tissue and hydrate until ready for next step
- 4. Block with 2% BSA (or blocking solution of choice) for 45 minutes
- 5. <u>1 Ab (or PBS for control) for 45 minutes</u>
- 6. Rinse 4X, 5 min. each with PBS
- 7. <u>2 Ab** (or PBS for control) for 1 hour</u>
- 8. Rinse 2X, 5 min. each with PBS
- 9. DAPI counterstain if applicable; 1 minute
- 10. Rinse 2X, 5 min. each with PBS
- 11. Mount coverslip; seal if necessary
- 12. Store covered (protected from light) at 4C. Temperature of storage may be dependent on fluorophore.
- * Each set of experiments should contain a negative control (the slide is subjected to the same procedure, but instead of adding antibodies all steps receive buffer) and a 2 Ab control (no primary is added, but 2 Ab is added.) If desired, a non-specific primary antibody control may also be performed.
- * All steps are performed at room temperature.
- ** Samples should be protected from light as much as possible after addition of 2 Ab. The easiest way to do this is to cover them with a box or other non-transparent object.